


FAST READ DIGITAL CLINICAL THERMOMETER

INSTRUCTIONS FOR USE


Please read carefully before using

KD-2050


Manufacturer :
K-Jump Health Co., Ltd
No.56 Wu Kung 5th Rd.,
Taipei, Taiwan

EC Representative :
PolyGreen Germany GmbH
Stahnsdorf, Germany


FEATURES

1. Fast oral and rectal temperature measurement in around 10 seconds with proper use.
2. ThermoTrend helps you better understand your health by showing the trend of last 3 temperatures measured.
2. Memory displays for the last temperature taken.
3. Very sensitive unit, best for quick oral or rectal measurement.
4. Easy to read digital LCD (liquid crystal display).
5. Compact, accurate and durable LSI (large scale integration) unit.
6. Temperature can be taken orally or rectally.
7. If the thermometer is inadvertently left on after temperature stabilization , it will automatically shut off in about 30 seconds.
8. Small , light weight unit . Universal use for the whole family. especially ideal for children.
9. ABS resin body presents no hazard compared to broken mercury /glass thermometers.
10. Low battery indicator : when the symbol "  "appears , the battery needs to be replaced .
11. This unit is waterproof

PARTS DESCRIPTION



SPECIFICATIONS


- Range :** 32.0~43.9°C(89.6~111.1°F)
Below 32.0°C(89.6°F) displays L°C(°F)
Above 43.9°C(111.1°F) displays H°C(°F)
- Resolution :** 0.1 °C(°F)
- Accuracy :** ±0.1°C (34.0°C to 42.0°C); ±0.2°C(other range)
±0.2°F(93.2°F to 107.6°F); ± 0.4°F(other range)
- Display :** Liquid crystal display 3 1/2 digits
- Battery : (included)** CR2032 (3V)
- Power Consumption :** 0.40 milliwatts in measurement mode
- Battery Life:** More than 160 hours of continuous operation
- Dimensions :** Approx. 117mm x 30mm x 16mm
- Weight :** Approx. 20grams including battery
- Beeper :** Approx. 10 seconds sound signal when peak temperature reached
- Memory :** Built-in memory enabling display of previously measured value
- Working Conditions :** Temperature: 10 ~ 40°C(50 ~ 104°F)
Relative Humidity: 15 ~ 95% non condensing
- Storage Conditions:** Temperature: -20~ 60°C(-4 ~ 140°F)
Relative Humidity: 15 ~ 95% non condensing
- Guarantee of Quality :** Certification ISO 13485
Comply with ASTM-E1112, EN12470-3, EN60601-1
- Product Classification :** Type BF  equipment




Temperature (°C)	Temperature (°F)	FeverTeller
T ≤35.9	T ≤96.7	Clear
36.0≤ T ≤36.5	96.8≤ T ≤97.7	Green
36.6≤ T ≤37.1	97.8≤ T ≤98.8	Light Green
37.2≤ T ≤37.4	98.9≤ T ≤99.4	Orange
37.5≤ T ≤37.9	99.5≤ T ≤100.3	Light Red
38.0≤ T ≤43.9	100.4≤ T ≤111.1	Red
T ≥ 44.0	T ≥ 111.2	Clear

PRECAUTION

1. For best results, do not use the fast read thermometer for axillary (underarm) measurement. Only use it for taking oral and rectal temperatures.
2. Do not boil the probe . Instead , clean the unit by wiping it with a dry cloth and disinfect the probe with ethyl (rubbing) alcohol .
3. Do not drop the thermometer or expose to heavy shock . The unit is not shockproof .
4. Do not bite the probe.
5. Do not store the unit in direct sunlight, or at a high temperature, or in high humidity or dust. Performance may be degraded.
6. Stop using the thermometer if it operates erratically or it displays malfunctions .
7. Keep out of the reach of unsupervised children .
8. Clean the thermometer probe before storing .
9. Do not attempt to disassemble the unit except to replace the battery .
10. If the thermometer has been stored at below-freezing temperatures, allow the thermometer to warm naturally to room temperature before using .
11. Performance of the device may be degraded if operated or stored outside stated temperatures and humidity ranges or the patient's temperature is below the ambient (room) temperature.

BATTERY REPLACEMENT



1. When the symbol "  " appears, the battery needs to be replaced.
2. Using a coin, unlock the battery cover by turning it counter-clockwise and carefully remove the battery with a non-metal instrument from the battery holder. Avoid using any sharp metal object in this operation.
3. Discard old battery carefully.
4. Replace the battery with a 3V, CR2032 lithium type or equivalent in the chamber with the negative (-) side down and the positive (+) side up.
5. Using a coin, replace and secure the battery cover by turning it clockwise.



Warning :

1. Discard old battery carefully , out of reach of young children. Swallowing the battery may be fatal. If the battery is swallowed, contact a hospital immediately to have it removed. Do not dispose of the battery in a fire. It may explode .
2. For battery disposition, please consider the national regulation.

HOW TO USE

1. Disinfect the probe with ethyl (rubbing) alcohol before using.
2. Depress the ON/OFF button. The display will read  or 
with FeverTeller flashing (Red→Orange→Green).
3. Release the power switch and the display will show L°C(L°F) with °C(°F) flashing.
4. Place the probe in the appropriate position (oral, axillary or rectal) .
5. Once the degree sign °C(°F) on the display stopped flashing , the correct temperature is indicated and FeverTeller will light up color according to the measured temperature. The temperature reading will not change after the °C(°F) stops flashing.
6. The unit will automatically turn off in 30 seconds (approx.). However, to prolong battery life, it is best to turn the thermometer off by pressing the ON/OFF button once the temperature has been noted.

Special Features:


Normal Temperature Alarm: An alarm will sound when the peak temperature has registered and is ready-to-read. The "Normal Temperature" alarm (four quick beeps followed by a pause) will sound repeatedly for approximately 10 seconds.

Last-Temperature Memory Recall: When the ON/OFF button is depressed to activate the thermometer, hold down for approximately 2 seconds. The display of 188.8 will change to show the last measured temperature. This reading is displayed as long as the button is held down. When the button is released, the thermometer will be in "ready-to-measure" mode and the recalled temperature will be erased.

ThermoTrend: It will show for 4 seconds (1) after last temperature recall (2) after measurement completed, with the last-temperature bar (right bar) blinking. The right bar shows the trend of last measurement, by comparing with 2nd last measurement, if last measurement is bigger (smaller) by 0.2°C(0.4°F) than 2nd last measurement, the bar will level up (down) one stage. The middle bar shows the trend of 2nd last measurement, by comparing with the 3rd last measurement. The left bar shows the trend of 3rd last measurement, by comparing with the 4th last measurement.

For example:

Temp Going up



Temp Going down



Temp at the same level



For Models with Fever Alarm : If the temperature registers higher than 37.5°C(99.5°F), the fever alarm (a quick beeping sound which occurs continuously for about 10 seconds) will sound when the peak temperature has registered. Also, the result will keep flashing until the thermometer is turned off.

For Switchable Centigrade and Fahrenheit Models : After turning on the thermometer, the Centigrade or Fahrenheit mode can be switched within the first 5 seconds by pressing the ON/OFF button.

NOTE : Always disinfect the thermometer before using.

NOTE : In Step 2, if the room temperature is greater than 32.0°C (89.6°F), it will be displayed rather than L°C(or L°F).

NOTE : In Step 3, you may put a sanitary plastic probe shield on the sensor tip to help prevent cross-infection and contamination. Use of the probe cover may result in a 0.1°C (0.2°F) difference from actual temperature. Dispose of any used probe shield and store the unit in a safe place.


RELIABLE METHODS TO MEASURE TEMPERATURE

THE TAKING OF BODY TEMPERATURES: Change in body temperature is one of the most important indicators of illness. It is important that accurate temperatures be taken. Even though the digital thermometer is used in similar manner to the glass thermometer, it is important that the simple instructions for use of this thermometer be followed closely and that the person taking the temperature understands body temperature.

BODY TEMPERATURE: A person's body temperature will usually be several tenths to one full degree lower in the morning than it is in the late afternoon. Therefore, if your temperature is 36.6°C(97.9°F) in the morning, it could be 36.9°C(98.4°F) or more in the late afternoon and still be normal. To determine what your normal body temperature is you should take your temperature both in the early morning and in the late afternoon over several days. Children have the same normal temperature range as do adults. A fever for both children and adults is typically defined as a temperature over 37.5°C(99.5°F) (oral reading). Exercise, excessive clothing, hot bath, hot weather, warm food and drink can cause mild elevations of up to 38°C (100.4°F) to 38.4°C(101.1°F) in oral temperature. If you suspect one of these to be the cause of an elevated reading, eliminate the possible cause and retake the temperature in a half hour.

ORAL USE

Place the probe well under the patient's tongue. Instruct the patient to keep their mouth closed while the thermometer is reading. A normal temperature by this method is usually considered between 36.0°C to 37.5°C(96.8°F to 99.5°F).



RECTAL USE

If rectal temperature is recommended by your doctor, you may use a probe shield. Insert thermometer into the probe shield. Lubricate the probe cover with a water soluble jelly for easier insertion . Do not use petroleum jelly . Insert tip of probe no more than 1/2 inch into rectum . STOP if you meet any resistance . The sensing unit is on the very tip of the probe and there is no need to insert the probe deep into the rectum . In general, the rectal temperature is around 0.5°C to 1.0°C(1.0 to 2.0°F) higher than oral temperature.

AXILLARY USE

Wipe armpit with a dry towel . Place probe in the patient's armpit and keep the patient's arm pressed firmly against the body. In general, the axillary temperature is around 0.5°C to 1.0°C(1.0 to 2.0°F) lower than oral temperature. Note that in order to achieve better axillary temperature measurement result , **a minimum measuring time of 3~5 minutes is recommended regardless of the beep sound.**

NOTE : Children should be attended by an adult during the entire temperature recording process . After use , the digital thermometer should be stored in a safe place, out of reach of young children. If you compare temperatures between two thermometers , it is possible to have two different readings even if both thermometers are within their accuracy tolerances. Differences in how you take each reading and the thermometer's tolerance range can account for nearly 0.1°C(0.2°F) variability.

CLEANING INSTRUCTIONS

To clean the thermometer, wash the tip with a solution of mild detergent and cool water. Disinfect the thermometer by wiping the sensor and lower stem with a cloth dipped in a household antiseptic solution such as rubbing alcohol.

LIMITED WARRANTY

This thermometer is guaranteed for one year from the date of purchase against manufacturer's defect under normal, household use. Following the directions closely will ensure years of dependable operation. If the thermometer does not function properly, first check the batteries (see battery replacement instructions). Replace it if necessary.

REV.091217

PN