


OMRON



Digital thermometer Model i-Temp


- Digitales Fieberthermometer
- Thermomètre digital
- Termometro digitale
- Termómetro digital
- Termómetro digital
- Digital termometer
- Digitalinen kuumemittari
- Digital termometer
- Digital febertermometer
- Digitale Koortsthermometer
- Ψυριακό Θερμόμετρο
- Dijital klinik termometresi
- Цифровой термометр

النحوت الرقمي

DE FR IT ES PT DK FI NO SE NL EL TR RU AR

IM-MC-670-E-03-030808
1664438-7C

A Good Sense of Health



Digital Thermometer English

| Indicator | Cause | Remedy |
|------------------------|----------------------------------------------------------------------------|----------------------------------------------------------------------------------|
| A Sensor tip | Displayed when the temperature of the sensing section is less than 32°C. | This is not an error so may be ignored. |
| L (low) | Displayed when the temperature of the sensing section is less than 32°C. | Check the sensing section and then take the temperature. |
| H (high) | Displayed when the temperature of the sensing section is higher than 42°C. | The thermometer is damaged. Please contact your nearest authorized OMRON dealer. |
| Eri - | The thermometer was stored at a temperature exceeding 40°C. | Allow the thermometer to adjust to room temperature before using it again. |
| RH - | The thermometer was stored at a temperature between 10°C and 40°C. | Take the thermometer to a temperature between 10°C and 40°C. |
| AL - | The thermometer was stored at a temperature below 10°C. | The thermometer cannot be used to measure the temperature in the rectal mode. |
| Air Temperature | Displayed when the surrounding temperature exceeds 32°C. | Check the sensing section and then take the temperature. |

5. Care and Maintenance

Keep the thermometer clean.
• Do not use any liquid to clean the thermometer. Do not allow water or other liquids to enter the unit.

Do not use any other method than cleaning and disinfecting. If you have any questions regarding specific cleaning methods, please refer to "Cleaning and Disinfecting" section of this manual.

Note: The thermometer cannot be used to measure the temperature in the rectal mode.

Turn off the thermometer if you hold it for 2 seconds or longer, the power will turn off. Turn off the power before returning it to its protective case.

Do not leave the thermometer in the storage case. Do not store the thermometer in the following locations:

- Locations with high heat and humidity or those that are exposed to direct sunlight.

- Areas where cleaning equipment, such as steam sterilizers, are used. Please refer to "Cleaning and Disinfecting" section of this manual.

Note: The thermometer was stored at a temperature exceeding 40°C.

Turn on the thermometer after leaving it at a temperature between 10°C and 40°C. Then take the temperature.

Do not leave the thermometer in a hot or cold location. Pharmaceutical storage areas or locations where chemicals are present.

Do not leave the thermometer in a car. - Wet locations.

Do not leave the thermometer in a car. - Locations with high heat and humidity or those that are exposed to direct sunlight.

Do not leave the thermometer in a car. - Areas where cleaning equipment, such as steam sterilizers, are used. Please refer to "Cleaning and Disinfecting" section of this manual.

Note: The thermometer was stored at a temperature exceeding 40°C.

Turn on the thermometer manually to conserve the battery. It will automatically switch off 1 minute after use when the power is turned off. Turn off the power before returning it to its protective case.

3. Battery replacement

Battery: Lithium Button Battery CR2016 (available commercially)

Product Name: MC-670-E

Power Supply: 3.0V DC, 1 CR2016 Lithium Button Battery

Power Consumption: 0.02 W

Battery Life: Approx. 1000 readings (if used once a day until the buzzer notification timing is reached)

Sensing Units: Thermistor

Measurement Method: Non-predictive

Measurement Site: Mouth (orally) or axilla (rectal)

Temperature Display: 3-digit, °C display in 1 degree increments

Measurement Accuracy: ±0.1°C (when measured at a standard room temperature of 23°C in a shaded area)

Measurement Range: 32.0 to 42.0°C


Usage Environment Temp and Humidity: Surrounding temperature 5 to 40°C, relative humidity 10% to 90%

Storage Environment Temp and Humidity: Surrounding temperature 20 to 60°C, relative humidity 10% to 95%

Height: Approx. 150 mm (6.3 inches)

External Dimensions: 34.0 mm (w) × 12.5 mm (h) × 5.0 mm (d)

Package Content: Test Battery & Lithium Button Battery CR2016, Storage Case, Instruction Sheet



Notes:

- This OMRON product is produced under the strict quality system of OMRON Healthcare Co.

- The specifications may be changed without prior notice.

- The accuracy is verified every two years by an OMRON representative.

Warranty:

This product is guaranteed by OMRON for 3 years after date of purchase. Any damage caused by incorrect handling shall not be covered by warranty.

Batteries and packaging are also excluded from the warranty. All other damage claims excluded.

A warranty claim must be submitted with the purchase receipt.



4. Read the Instructions in this manual carefully.

This device falls into the provisions of the EC directive 93/42/EEC (Medical Device Directive) and the European Standard EN12470:2000, Clinical thermometers - Part 3: Performance of compact electrical thermometers (non-predictive and predictive) with maximum device.

5. Use of the thermometer

1. Remove the thermometer from the storage case and press the ON/OFF button. (Fig. 5)

The thermometer will start beeping when it is turned on for a few seconds.

2. Turn the thermometer to the mouth (orally) or armpit (rectal).

3. Press the ON/OFF button again to stop the beeping.

4. Open the cover of the compartment for the batteries.

5. Insert the batteries.

6. Close the cover of the compartment for the batteries.

7. Turn the thermometer to the mouth (orally) or armpit (rectal).

8. Press the ON/OFF button again to start the beeping.

9. Turn the thermometer to the mouth (orally) or armpit (rectal).

10. Press the ON/OFF button again to stop the beeping.

11. Turn the thermometer to the mouth (orally) or armpit (rectal).

12. Press the ON/OFF button again to start the beeping.

13. Turn the thermometer to the mouth (orally) or armpit (rectal).

14. Press the ON/OFF button again to stop the beeping.

15. Turn the thermometer to the mouth (orally) or armpit (rectal).

16. Press the ON/OFF button again to start the beeping.

17. Turn the thermometer to the mouth (orally) or armpit (rectal).

18. Press the ON/OFF button again to stop the beeping.

19. Turn the thermometer to the mouth (orally) or armpit (rectal).

20. Press the ON/OFF button again to start the beeping.

21. Turn the thermometer to the mouth (orally) or armpit (rectal).

22. Press the ON/OFF button again to stop the beeping.

23. Turn the thermometer to the mouth (orally) or armpit (rectal).

24. Press the ON/OFF button again to start the beeping.

25. Turn the thermometer to the mouth (orally) or armpit (rectal).

26. Press the ON/OFF button again to stop the beeping.

27. Turn the thermometer to the mouth (orally) or armpit (rectal).

28. Press the ON/OFF button again to start the beeping.

29. Turn the thermometer to the mouth (orally) or armpit (rectal).

30. Press the ON/OFF button again to stop the beeping.

31. Turn the thermometer to the mouth (orally) or armpit (rectal).

32. Press the ON/OFF button again to start the beeping.

33. Turn the thermometer to the mouth (orally) or armpit (rectal).

34. Press the ON/OFF button again to stop the beeping.

35. Turn the thermometer to the mouth (orally) or armpit (rectal).

36. Press the ON/OFF button again to start the beeping.

37. Turn the thermometer to the mouth (orally) or armpit (rectal).

38. Press the ON/OFF button again to stop the beeping.

39. Turn the thermometer to the mouth (orally) or armpit (rectal).

40. Press the ON/OFF button again to start the beeping.

41. Turn the thermometer to the mouth (orally) or armpit (rectal).

42. Press the ON/OFF button again to stop the beeping.

43. Turn the thermometer to the mouth (orally) or armpit (rectal).

44. Press the ON/OFF button again to start the beeping.

45. Turn the thermometer to the mouth (orally) or armpit (rectal).

46. Press the ON/OFF button again to stop the beeping.

47. Turn the thermometer to the mouth (orally) or armpit (rectal).

48. Press the ON/OFF button again to start the beeping.

49. Turn the thermometer to the mouth (orally) or armpit (rectal).

50. Press the ON/OFF button again to stop the beeping.

51. Turn the thermometer to the mouth (orally) or armpit (rectal).

52. Press the ON/OFF button again to start the beeping.

53. Turn the thermometer to the mouth (orally) or armpit (rectal).

54. Press the ON/OFF button again to stop the beeping.

55. Turn the thermometer to the mouth (orally) or armpit (rectal).

56. Press the ON/OFF button again to start the beeping.

57. Turn the thermometer to the mouth (orally) or armpit (rectal).

58. Press the ON/OFF button again to stop the beeping.

59. Turn the thermometer to the mouth (orally) or armpit (rectal).

60. Press the ON/OFF button again to start the beeping.

61. Turn the thermometer to the mouth (orally) or armpit (rectal).

62. Press the ON/OFF button again to stop the beeping.

63. Turn the thermometer to the mouth (orally) or armpit (rectal).

64. Press the ON/OFF button again to start the beeping.

65. Turn the thermometer to the mouth (orally) or armpit (rectal).

66. Press the ON/OFF button again to stop the beeping.

67. Turn the thermometer to the mouth (orally) or armpit (rectal).

68. Press the ON/OFF button again to start the beeping.

69. Turn the thermometer to the mouth (orally) or armpit (rectal).

70. Press the ON/OFF button again to stop the beeping.

71. Turn the thermometer to the mouth (orally) or armpit (rectal).

72. Press the ON/OFF button again to start the beeping.

73. Turn the thermometer to the mouth (orally) or armpit (rectal).

74. Press the ON/OFF button again to stop the beeping.

