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# FlexTemp II English

**1. IMPORTANT SAFETY INSTRUCTIONS**

- The thermometer is suitable exclusively for measurement of the body temperature!
- The minimum measurement time until the signaling tone (bip) must be maintained without exception!
- The thermometer contains small parts (battery, etc.) that can be swallowed by children. Keep the thermometer out of the reach of children.
- Damage may occur if the thermometer is dropped or struck against a hard surface.
- Avoid ambient temperatures above 60°C. NEVER boil the thermometer!
- Use only the commercial disinfectants listed in the section «cleaning and disinfecting» to clean the thermometer. The thermometer must be intact when immersed in liquid disinfectant.
- **Warning:** The thermometer contains a long-life battery. The (■) symbol in the right bottom corner of the display indicates the battery needs to be changed. Open the battery compartment (■) by pulling the colored top firmly. Gently pull out the circuit board with battery compartment approximately 1 cm. Remove the battery by pushing a thin plastic or wooden tip (such as toothpick) through the small gap at the side of the battery holder. Push further to remove the battery from the holder. Replace the used battery with a new one paying attention to the "+" symbol on the battery which must be positioned upwards. For disposal, please comply with local regulations for special waste.
- The manufacturer recommends that the accuracy is verified by an authorized laboratory every 2 years. Please contact our Customer Services (0870 500 2771 0845) for advice (for UK only).

**2. SWITCHING ON**  
 Press the on/off button (■). A short beep-signal sounds, the device is on. All symbols show on the LCD display (■) for a few seconds as a means of checking that the thermometer is functioning correctly. If available the last correctly measured temperature is displayed indicated by the "M" symbol. After this the device is ready for measurement indicated by a flashing "°C" symbol and a "Lo" indication at an ambient temperature of less than 32°C.

**3. HOW TO USE THE THERMOMETER**  
 Position the sensor tip (■) in the preferred place of measurement - mouth, rectum or armpit. The signaling tone (beep-beep-beep) sounds 4 times and the «°C» no longer flashes, this means that the measurement is finished and that the thermometer is ready for reading.

**In the mouth (oral)**  
 Position the thermometer under the tongue. The measuring sensor must be in good contact with the tissue. Close your mouth and breathe evenly through the nose. Do not swallow. Wait until the measurement has been influenced by inhaled/exhaled air. Approx. measuring time: 1 minute. Normal temperature between 35.7 °C and 37.3 °C.

**In the anus (rectal)**  
 This is the most reliable measuring method and is especially appropriate for infants and small children. Carefully insert the tip of the thermometer 2 to 3 cm into the anal aperture. Approx. measuring time: 1 minute. Normal temperature between 36.2 °C and 37.7 °C.

**In the armpit (axillary)**  
 From a medical viewpoint, this method will always provide inaccurate readings and should not be used if precise measurements are required. Approx. measuring time: 1.5 minutes. Normal temperature between 35.2 °C and 36.7 °C.

**Warning!** To prolong the life of the battery, turn it off after use by pressing the operating button briefly. Otherwise the thermometer will automatically turn off after about 10 minutes.

| Manufacturer                 | Subsidiary             |
|------------------------------|------------------------|
| OMRON HEALTHCARE EUROPE B.V. | OMRON Santé France SAS |
| OMRON HEALTHCARE UK LTD.     | OMRON Santé France SAS |

| Manufacturer                 | Subsidiary             |
|------------------------------|------------------------|
| OMRON HEALTHCARE EUROPE B.V. | OMRON Santé France SAS |
| OMRON HEALTHCARE UK LTD.     | OMRON Santé France SAS |

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